WGS of Cholera: Building Global Capacity for Molecular Subtyping of Enteric Bacteria through PulseNet International

15 April 2019 – Annecy, France

Maryann Turnsek and Michele Parsons Centers for Disease Control and Prevention, Atlanta USA

Objectives

- Provide an overview of approaches and activities in PulseNet International to sequence enteric bacteria, including Vibrio cholerae
- Discuss challenges for implementation and data sharing regarding WGS

Overview of approach for sequencing Vibrio cholerae

- Applications for sequence data
 - Outbreaks
 - Surveillance
 - Research/special studies
- Preferred methods for WGS
 - Illumina chemistry and platforms
 - BioNumerics 7.6 or higher for analysis
- Managing sequence data
 - Raw reads maintained by CDC
 - Data uploaded in real time to Sequence Read Archive (SRA) at NCBI
 - Fastq files, limited metadata
 - Vibrio bioproject identifier: PRJNA266293









BioNumerics version 7.6





- Assembly
- Trimming, quality checks
- Average Nucleotide Identity (ANI)

- 'O' group determination
- Serotype, biotype
- Virulence profiling
- Resistance prediction
- cgMLST/wgMLST for surveillance

✓ Vibrio allele database is being validated and will be finalized in 2020.

Why use MLST to analyze and compare bacterial genomes?

- Stable nomenclature
- Computationally easier and faster than other methods
- No need for a reference strain





Why use BioNumerics Software for WGS data?

- Standardized way to manage and analyze sequence data
- 'One shot' reference characterization within customized organism-specific databases
- Does not require advanced bioinformatics skills





Customized for NCBI Submissions

NCBI submission settings						
BioProject and organisation Provide general information for all submissions: what BioProject do they belong to, and who is submitting.						
BioProject accession:	PRJNA218110					
Organisation name:	EDLB-CDC					
SPUID namespace:	EDLB-CDC					
Туре:	consortium					
Role:	owner 👻					
Contact first name:	EDLB					
Contact last name:	CDC					
Contact e-mail:	pulsenet@cdc.gov					
FTP upload directory:	submit/Test					
	< Back Next > Cancel					

BioProject accessions: Campylobacter: PRJNA239251 E. coli and Shigella: PRJNA218110 Listeria: PRJNA212117 Salmonella: PRJNA230403 Vibrio: PRJNA266293

PulseNet International – Implementing WGS for Foodborne Disease Surveillance

Why are PulseNet International lab moving towards WGS?

- WGS is more precise and provides more information than current "gold standard" subtyping techniques, including pulsed-field gel electrophoresis (PFGE) and multi-locus variable number tandem repeat analysis (MLVA)
- Data can be shared across all laboratories for routine surveillance, outbreak identification, antimicrobial resistance prediction and reference characterization

https://www.cdc.gov/pulsenet/participants/international/wgs-vision.html

PulseNet International



PulseNet International WGS Trainings

• Sequencing on the MiSeq Platform – Thailand, 2019





• Analysis methods – Kenya, 2019



Summary of WGS and Analysis Trainings Lead by CDC

Date	WGS wet lab	BioNumerics 7.6	Training site	Countries	#
November 2015		Х	PNLAC meeting, La Paz Mexico	Chile, Paraguay, Costa Rica, Mexico, Peru and Columbia, Argentina, Ecuador	25
March 2016	Х	Х	CDC-Atlanta	Chile	1
November 2016		Х	CDC-Atlanta	Costa Rica	3
March 2017	Х		InDRE	Mexico, Cuba	10
March 2017		Х	InDRE	Mexico	5
June 2017	Х	X	CDC-Atlanta	Kenya	1
July-August 2017	Х	Х	Chile Institute of Public Health	Chile, Paraguay, Brazil, Peru and Argentina	9
November 2017	Х	Х	CDC-Atlanta	Argentina	2
June 2018	Х	Х	CDC-Atlanta	Lebanon and Israel	2
July 2018	Х	Х	KEMRI	Kenya	18
February 2019	Х	Х	CDC-Thailand	Thailand	15

Challenges Implementing WGS in Countries with Highest Cholera Burden

- Availability of reagents limited ordering options (usually only regional distributors), procurement delays common
- Instrument care sequencing platforms are sensitive, maintenance contracts are highly recommend, but expensive and require highly trained technicians
- IT infrastructure improvements needed to handle / store volume of sequence data generated, ensure security of data
- Availability of analytical tools limited "user friendly" and affordable options
- Internet speeds slow internet speeds limit analysis, data transfer
- High cost depending on Illumina platform and chemistry (capacity ranges from 2 40 genomes), \$1 000 \$4 000 USD / run, \$100 \$500 USD / sample (not including cost to analyze or store data)

Data sharing

- At trainings, we encourage standard policy for data sharing
 ✓ Real time release to SRA in NCBI with limited metadata
- For requests, we share interpretation (e.g. phylogenetic tree)
 ✓ Important to see what the data looks like
- From our experience, there is a great deal of mis-trust of genetic data being in the public realm
 - ✓ Lead by example

CDC Laboratory Contacts for GTFCC PulseNet International Contacts

- Michele Parsons, PhD
 - ypi5@cdc.gov
- Maryann Turnsek, BS
 - hud4@cdc.gov
- Cheryl Tarr, PhD
 - crt6@cdc.gov

- Jeny Concepción-Acevedo, PhD
 - Coordinator, PulseNet International
 - ypi5@cdc.gov
- Peter Gerner-Smidt, MD PhD
 - Chair, PulseNet International Steering Committee
 - plg5@cdc.gov

Acknowledgements:

- Enteric Diseases Laboratory Branch (EDLB) in the Division of Foodborne, Waterborne and Environmental Diseases, National Center for Emerging and Zoonotic Infectious Diseases, CDC
- Epidemiology, Informatics, Surveillance and Laboratory Branch (EISLB) in the Division of Global Health Protection, Center for Global Health, CDC



The findings and conclusions in this report are those of the authors and do not necessarily represent the official position of the Centers for Disease Control and Prevention.

